*plastlabor*

10

Microbiologics

**FOOD SAFETY  
CONTROLS**

**FOOD SAFETY QC MICROORGANISMS**

**Microbiologics** 

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## Recommended Culture Methods for Microorganisms

### Selection of Growth Requirements

1. Primary growth on a nonselective agar medium is preferred. Primary growth in a fluid medium should only occur in special instances or when recommended. Because of the manipulations required during hydration, it is difficult to obtain purity of a lyophilized strain in a fluid medium. A contaminant may completely overgrow and obscure the presence of the lyophilized strain.
2. The following information lists which method should be used to grow the various microorganism species. Descriptions of methods follow the microorganism list.

Microorganism	Method	Notes
<i>Acetobacter</i> species	Method 36	
<i>Achromobacter</i> species	Method 1	
<i>Acinetobacter</i> species	Method 1	
<i>Actinobacillus</i> species	Method 3	
<i>Actinomyces</i> species	Method 4	
<i>Aerococcus</i> species	Method 1	
<i>Aeromonas</i> species	Method 2	Exceptions are <i>Aeromonas hydrophila</i> , Microbiologics 0870 and <i>Aeromonas salmonicida</i> .
<i>Aeromonas hydrophila</i> , Microbiologics 0870	Method 31	
<i>Aeromonas salmonicida</i>	Method 32	
<i>Aggregatibacter</i> species	Method 3	
<i>Alcaligenes</i> species	Method 1	
<i>Allicyclobacillus</i> species	Method 12	An exception is <i>Allicyclobacillus acidoterrestris</i> .
<i>Allicyclobacillus acidoterrestris</i>	Method 45	
<i>Allofococcus</i> species	Method 2	
<i>Atemaria</i> species	Method 5	
<i>Aspergillus</i> species	Method 5	
<i>Anaerobaculum</i> species	Method 1	
<i>Aquaspirillum</i> species	Method 20	
<i>Arcanobacterium</i> species	Method 34	
<i>Arthrobacter</i> species	Method 21	
<i>Aspergillus</i> species	Method 5	An exception is <i>Aspergillus flavus</i> .
<i>Aspergillus flavus</i>	Method 46	
<i>Aureobasidium</i> species	Method 5	
<i>Bacillus</i> species	Method 49	
<i>Bacteroides</i> species	Method 4	
<i>Bifidobacterium</i> species	Method 4	An exception is <i>Bifidobacterium animalis</i> subsp. <i>animalis</i> .
<i>Bifidobacterium animalis</i> subsp. <i>animalis</i>	Method 39	
<i>Bordetella bronchiseptica</i>	Method 15	

3. The following information lists methods for growing microorganisms. When possible, more than one type of agar medium per method is listed.

#### Method 1

Media	Tryptic Soy Agar (Soybean Casein Digest Agar), Nonselective Sheep Blood Agar, Standard Methods Agar (Plate Count Agar) or Nutrient Agar
Temperature	35°C
Atmosphere	Aerobic
Growth Time	24 to 48 hours

#### Method 2

Media	Nonselective Sheep Blood Agar
Temperature	35°C
Atmosphere	Aerobic
Growth Time	24 to 72 hours

#### Method 3

Media	Chocolate Agar
Temperature	35°C
Atmosphere	5 to 7% Carbon Dioxide
Growth Time	24 to 48 hours

#### Method 4

Media	Anaerobic Blood Agar
Temperature	35°C
Atmosphere	Anaerobic
Growth Time	48 to 72 hours

Note: Some obligate anaerobes may require 5-7 days to demonstrate sufficient growth.

#### Method 5

Media	Sabouraud Dextrose Emmons Agar
Temperature	25°C
Atmosphere	Aerobic
Growth Time	2 to 7 days

Note: Nonselective Sheep Blood Agar is an appropriate alternative. Nutrient Agar, Tryptic Soy Agar, Potato Dextrose Agar, and Standard Methods Agar (Plate Count Agar) are appropriate alternatives together with an additional period (24 hours) of incubation.

#### Method 6

Media	Chocolate Agar
Temperature	35°C
Atmosphere	Microaerophilic
Growth Time	48 to 72 hours

Note: Do not open the inoculated agar medium petri plate for the first 48 hours.

## Microorganism Maintenance Plan for Laboratories Following ISO 11133

KWIK-STIK™ and LYFO DISK™ products are reference stock cultures derived from reference cultures and preserved by a process called lyophilization. According to ISO 11133, 5.3.3, reference stock cultures are to be “properly maintained and handled in a manner that minimizes the opportunity for cross-contamination, mutation, or alteration of typical characteristics.”

This technical information bulletin provides a month-long maintenance plan for preserving the viability, purity, and genotypic and phenotypic characteristics of a microorganism strain grown from KWIK-STIK or LYFO DISK. The main plan is based on the plan described in Chart B.2 of ISO 11133<sup>1</sup>.

### Getting Started

- After hydration, plate the KWIK-STIK or LYFO DISK on a non-selective agar, and incubate under conditions suitable for the microorganism strain.
- If you are unsure what agar to use, there are three places to find culture methods:
  - ISO 11133, Tables E and F.
  - Microbiologics website, [www.microbiologics.com](http://www.microbiologics.com). Recommended culture methods are listed on each strain's webpage.
  - The Technical Information Bulletin, *Recommended Culture Methods*, which is available at [www.microbiologics.com](http://www.microbiologics.com)

### Storage of Reconstituted Microorganisms

- Most quality control microorganisms can be maintained on nonselective agar plates or slants for up to four weeks at room temperature or in the refrigerator.
- Microbiologics has found the following storage conditions to be favorable for maintenance:

Category of Microorganism	Storage Conditions
Aerobic Bacteria	Store at 2-8°C. A few species of <i>Bacillus</i> remain viable for a longer period when stored at room temperature.
Yeast and Fungi	Store at room temperature.
Anaerobes	Store in anaerobic conditions at room temperature.
<i>Campylobacter</i> species	Store on chocolate agar at 35°C in microaerophilic conditions.

### Five Reasons Not to Freeze


- It is not necessary! Microbiologics has done the hard work of confirming the identity and purity of a strain and preserving it for you.
- Freezing adds a passage. ISO 11133, 5.3.1, says, “The number of transfers of test organisms should be documented to prevent excessive sub-culturing that increases the risk of phenotypic alteration. One passage is defined as the transfer from a viable culture to a fresh medium with growth of the microorganisms. Any form of sub-culturing is considered to be a form of transfer/passage.”
- A -70°C freezer is needed. ISO 11133, 5.3.3, says, “Reference stocks should be stored in multiple portions, usually either deep-frozen, e.g. below -70°C, or lyophilized. At a higher temperature, duration of viability might be reduced and genetic modification might occur.”
- Microbiologics cannot guarantee the typical characteristics of a microorganism if it is frozen. We have not validated the stated characteristics remain the same after the microorganism has been frozen and thawed.
- ISO 11133, 5.3.2, says, “If test microorganisms are obtained from reference collections or commercial suppliers holding ISO 9001 certification or other appropriate certification and maintained in their original containers, the manufacturer's directions for their cultivation and use shall be followed.”

REFERENCES  
ISO 11133:2014 Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media

### Tips for Best Performance

- Do not test the original pellet plate for phenotypic characteristics. This is the plate on which the lyophilized pellet was started. The organisms growing on this plate are not fully resuscitated.
- Select isolated colonies for a test. Do not test colonies from a contaminated plate.
- When possible, use microorganisms that are not more than 24 hours old for biochemical tests.
- Microorganisms stored at 4°C should not be used for certain tests (e.g. MALDI-TOF). Consult the manufacturer's instructions.
- After the fourth week, dispose of plates and start the process over with a new lyophilized pellet.
- A microorganism may be used beyond expiry date if (1) the lyophilized pellet is grown before expiry date and (2) the microorganism is not used beyond week four of the maintenance program.

## Microorganism Maintenance Plan for Strains Used in ISO 11133

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## BIG 6

## BIG 7

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- Available in quantitative and qualitative formats<sup>1</sup>
- Ensures accurate test results
- Supports compliance with USDA/FDA regulations
- Available internationally<sup>1</sup>

STEC Strain	KWIK-STIK™	LYFO DISK™	Epower™	UN-BioTAG™
01100 <i>E. coli</i> serotype O26:H11 derived from CDC 03-3014	•	•	•	
01098 <i>E. coli</i> serotype O45:H2 derived from CDC 00-3039	•	•	•	
01101 <i>E. coli</i> serotype O103:H11 derived from CDC 06-3008	•	•	•	
01102 <i>E. coli</i> serotype O111:H8 derived from CDC 2010C-3114	•	•	•	
01099 <i>E. coli</i> serotype O121:H19 derived from CDC 02-3211	•	•		
01097 <i>E. coli</i> serotype O145:NM derived from CDC 99-3311	•	•	•	
0617 <i>E. coli</i> serotype O157:H7 derived from ATCC® 35150™*	•	•	•	
01104 <i>E. coli</i> serotype O104:H4 derived from ATCC® BAA-2326™*	•	•	•	
01227 <i>E. coli</i> serotype O157:H7 EC43 derived from FDA ESC1177				•

QC Sets and Panels are also available: 5218 Big 6 QC Set and 5219 Big 7 STEC QC Set

<sup>1</sup>Epower™ and UN-BioTAG™ are not intended for clinical use. \*Restrictions may apply by region.

ATCC Licensed  
\*Look for the ATCC Licensed Derivative® symbol for products derived from ATCC® cultures. The ATCC Licensed Derivative symbol, the ATCC Licensed Derivative Word mark, and the ATCC Celling Mark are trademarks of ATCC - Microbiology, Inc. In licensed to use these trademarks and to sell products derived from ATCC® cultures.

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**UV-BioTAG**

Catalog # / Format	Strain	Characterization
01227UV-S / Swab Kit 01227UV-V / Vial Kit	<i>Escherichia coli</i> (O157:H7) EC43 (STEC) derived from FDA EXC 11177 (STEC)	Chloramphenicol (Cm10 <sup>6</sup> ) resistant
01223UV-S / Swab Kit 01223UV-V / Vial Kit	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Typhimurium SAL54 derived from FDA SAL 5694 / FDA SAL 5695	H <sub>2</sub> S+ Chloramphenicol (Cm10 <sup>6</sup> ) resistant
01277UV-S / Swab Kit 01277UV-V / Vial Kit	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Typhimurium Sal46 derived from FDA SAL5693	H <sub>2</sub> S + Chloramphenicol (Cm10 <sup>6</sup> ) resistant
01226UV-S / Swab Kit 01226UV-V / Vial Kit	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Settenberg SAL59 derived from FDA SAL 5697	Lactose - H <sub>2</sub> S - Chloramphenicol (Cm10 <sup>6</sup> ) resistant
01278UV-S / Swab Kit 01278UV-V / Vial Kit	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Gammaria Sa67 derived from FDA SAL5695	H <sub>2</sub> S + Chloramphenicol (Cm10 <sup>6</sup> ) resistant
01230UV-S / Swab Kit 01230UV-V / Vial Kit	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Mbandaka Sa58 derived from FDA SAL5656	Sucrose + H <sub>2</sub> S + Chloramphenicol (Cm10 <sup>6</sup> ) resistant
01228UV-S / Swab Kit 01228UV-V / Vial Kit	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Minnesota Sa50 derived from FDA SAL5658	Lactose - H <sub>2</sub> S + Chloramphenicol (Cm10 <sup>6</sup> ) resistant
01224UV-S / Swab Kit 01224UV-V / Vial Kit	<i>Shigella flexneri</i> SF64 derived from FDA SH6552 / FDA SH6553	Kanamycin (Km25 <sup>6</sup> ) resistant Chloramphenicol (Cm10 <sup>6</sup> ) resistant
01276UV-S / Swab Kit 01276UV-V / Vial Kit	<i>Shigella sonnei</i> S561 derived from FDA SH6551	Chloramphenicol (Cm10 <sup>6</sup> ) resistant
01249UV-S / Swab Kit 01249UV-V / Vial Kit	<i>Listeria monocytogenes</i> (1/2a) derived from FDA LS808	Serotype 1/2a Erythromycin resistant
01248UV-S / Swab Kit 01248UV-V / Vial Kit	<i>Listeria monocytogenes</i> (1/2b) derived from FDA LS810	Serotype 1/2b Erythromycin resistant Rifampicin resistant
01250UV-S / Swab Kit 01250UV-V / Vial Kit	<i>Listeria monocytogenes</i> (4b) derived from FDA LS806	Serotype 4b Erythromycin resistant Streptomycin resistant
01271UV-S / Swab Kit 01271UV-V / Vial Kit	<i>Listeria innocua</i> (5a) modified from NCTC 11288	Serotype 5a Erythromycin resistant

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# Como utilizar?

**1**

Allow the unopened UV-BioTAG swab pouch to equilibrate to room temperature. Tear open pouch at notch and remove the UV-BioTAG swab unit.

**2**

Break red snap valve at the top of the UV-BioTAG swab to release the hydrating fluid.

**3**


Squeeze the bulb at the top of the UV-BioTAG swab to rehydrate the pellet.

**4**

Using a pinching action on the bottom portion of the unit, crush the pellet in the fluid until the pellet suspension is homogenous.

**5**

Inoculate a primary culture plate(s) by gently rolling the swab over one-third of the plate.

**6**

Continue according to laboratory protocol.

**1**

Remove the plastic container containing the vials of pellets from refrigerated storage. Remove the vials to be used; immediately place the plastic container containing the remaining vials back into refrigerated storage to maintain product stability.

**2**

Aseptically remove 1 pellet with sterile forceps from the vial. Do not remove desiccant.

**3**

Place the pellet in 0.5 ml of sterile fluid (water, saline, Tryptic Soy Broth, or Brain Heart Infusion Broth).

**4**

Crush the pellet with a sterile swab until the suspension is homogenous.

**5**

Inoculate a primary culture plate(s) by gently rolling the swab over one-third of the plate.

**6**

Continue according to laboratory protocol.

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Certificate of Analysis: Lyophilized Microorganism Specification and Performance Upon Release

<b>Specifications</b> Microorganism Name: <i>Listeria monocytogenes</i> (1/2a) Catalog Number: 01249 Lot Number: 1249-03** Reference Number: FDA LS808 Purity: Pure Passage from Reference: 3	<b>Expiration Date:</b> 2020/7/31 <b>Release Information:</b> Quality Control Technologist: Mary L. Bowman Release Date: 2018/8/30
---	---

Performance	
<b>Macroscopic Features:</b> Small, circular, pale white, translucent.	<b>Medium:</b> SBAP
<b>Microscopic Features:</b> Gram positive rods.	<b>Method:</b> Gram Stain (1)
<b>ID System:</b> MALDI-TOF (1) See attached ID System results document.	<b>Other Features/ Challenges: Results</b> (1) Catalase (3% Hydrogen Peroxide): positive (1) Purple Broth w/Rhamnose: positive (1) Purple Broth w/Rhamnose: negative <b>Longwave UV fluorescence: fluorescent to five passages from pellet best visualized using blue light excitation 475 - 495 nm light and blue light barrier lenses.</b>

Longwave UV fluorescence: fluorescent to five passages from pellet best visualized using blue light excitation 475 - 495 nm light and blue light barrier lenses.

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
Let's Get Growing


## Microbiologics Dilutions Guide

This document outlines how to perform dilutions when using Microbiologics products. **Dilution is the process of making a solution weaker or less concentrated.** In microbiology, serial dilutions (log dilutions) are used to decrease a bacterial concentration to a required concentration for a specific test method, or to a concentration which is easier to count when plated to an agar plate. This document was created to provide a better understanding of dilutions and should be used as a guideline, not a replacement for laboratory procedures.

### Log Dilutions

A log dilution is a tenfold dilution, meaning the concentration is decreased by a multiple of ten. To complete a tenfold dilution, the ratio must be 1:10. The 1 represents the amount of sample added. The 10 represents the total size of the final sample. For example, a sample size of 1 ml is added to 9 ml of diluent to equal a total of 10 ml.

**Example:** 1:10 dilution - if the concentration is 1,000 CFU, a one log dilution will drop the concentration to 100 CFU.

### Decimal Numbers vs Scientific Notation

Decimal numbers can be converted to scientific notations by moving the decimal place the same number of places as the exponential number.


$$100.0 = 1.0 \times 10^2 \quad 1.0 \overbrace{00}^2 = 100$$

$$350.0 = 3.5 \times 10^2 \quad 3.5 \overbrace{00}^2 = 350$$

Decimal Number	1	10	100	1,000	10,000	100,000	1,000,000
Scientific Notation	$10^0$	$10^1$	$10^2$	$10^3$	$10^4$	$10^5$	$10^6$

### Multiple Dilutions

Multiple dilutions are required to decrease the sample concentration by multiple logs. If the concentration is 35,000 CFU/ml ( $10^4$ ), and 35 CFU/ml is the target concentration, the following serial dilutions can be performed.


### Larger Dilutions

Decreasing the concentration using fewer dilutions is possible with the use of large volume dilutions. This can be done by performing a 1:100 dilution instead of 1:10. An example of this can be observed in the Epower™ instructions for membrane filtration using Microbiologics E3 Epower™ product.

The E3 Epower™ product provides  $10^3$  CFU per pellet which equates to 1 pellet in 1 ml equaling  $10^3$  CFU/ml. Placing a  $10^3$  Epower™ pellet in 10 ml will drop the concentration to  $10^2$  - this is a one log dilution.

A 1:100 dilution can be created by placing 1 pellet in 99 ml as instructed in the membrane filtration instructions. This will drop the concentration two logs from  $10^3$  to  $10^1$  CFU/ml.

### EZ-CFU™ Dilution Example

When EZ-CFU™ is used according to directions, the following dilutions are conducted to reach a desired concentration of 10-100 CFU/0.1 ml.

- Two pellets are placed in 2 ml of hydrating fluid = 1,000 - 10,000 CFU/ml.
- 1:10 dilution is performed by placing 1 ml of the re-hydrated pellet solution into 9 ml of buffer = 100 - 1,000 CFU/ml.
- 0.1 ml of the organism suspension plated to an agar = 10 - 100 CFU per 0.1 ml.

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Pellet Concentration	Examples of Concentration (CFU/ml) in Specified Hydrating Fluid Volume		
	1 ml	10 ml	100 ml
10,000 – 99,000	10,000 – 99,999	1,000 – 9,999	100 - 999

- 1

Remove the vial of pellets from refrigerated storage. Allow the unopened vial to equilibrate to room temperature (about 30 minutes).
- 2

Prior to use, warm hydrating and dilution fluids to 34°C-38°C. Sterile pH 7.2 phosphate buffer is recommended for hydration of the lyophilized preparation.
- 3

With a sterile forceps, transfer the **EZ-SPORE™ microorganism** pellet to the hydrating fluid. Do not remove the desiccant from vial. Immediately stopper and recap vial and return to 2°C-8°C.
- 4

Place the microorganism suspension into a 34°C-38°C incubator for 30 minutes to ensure complete hydration.
- 5

Immediately following incubation, mix hydrated material until a homogeneous suspension is achieved.
- 6

Proceed with the challenge according to laboratory protocol.
- 7

The challenge must be completed within 30 minutes of the hydration process to avoid a change in the challenge suspension concentration.

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**Epower CRM (ISO 17034)**

Certificate of Analysis: Lyophilized Microorganism Specification and Performance Upon Release

<p><b>Specifications</b>  <b>Microorganism Name:</b> Clostridium perfringens  <b>Catalog Number:</b> 0318(CRM)(6)  <b>Lot Number:</b> 318-236**  <b>Reference Number:</b> ATCC® 13124™*  <b>Passage from Reference:</b> 1</p>	<p><b>Expiration Date:</b> 2022/3/31  <b>Release Information:</b>  <b>Quality Control Technologist:</b> Kavitha Gobalan  <b>Release Date:</b> 2020/5/14</p>
---	---

Certificate of Analysis: Lyophilized Microorganism Specification and Performance Upon Release

<p><b>Specifications</b>  <b>Microorganism Name:</b> Escherichia coli  <b>Catalog Number:</b> 0495(CRM)  <b>Lot Number:</b> 495-123**  <b>Reference Number:</b> ATCC® 35218™*  <b>Passage from Reference:</b> 1</p>	<p><b>Expiration Date:</b> 2022/3/31  <b>Certified Value:</b> 6.0E+03 CFU per pellet  <b>(2) Expanded Uncertainty:</b>            95% Confidence Interval: 2.8E+03 to 1.3E+04 CFU  <b>(3) Standard Deviation:</b> 4.7E+02  <b>Release Information:</b>  <b>Quality Control Technologist:</b> Christine Condon  <b>Release Date:</b> 2020/4/30</p>
---	---

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