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Article - Human and Animal Health

Detection of Drug Resistant Mycobacterium Tuberculosis Strains Using Kit SIRE Nitratase®: a Multicenter Study

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Miguel Viveiros¹⁴
<https://orcid.org/0000-0001-9676-6251>

Afrânio Lineu Kritski¹
<https://orcid.org/0000-0002-5900-6007>

¹Federal University of Minas Gerais, Faculty of Medicine, Mycobacteria Research Laboratory, Belo Horizonte, Minas Gerais, Brazil; ²Federal University of Minas Gerais, Faculty of Pharmacy, Department of Social Pharmacy, Belo Horizonte, Minas Gerais, Brazil; ³Federal University of Minas Gerais, Veterinary School, Department of Preventive Veterinary Medicine, Belo Horizonte, Minas Gerais, Brazil; ⁴Federal University of Rio Grande, Faculty of Medicine, Laboratory of Mycobacteria, Rio Grande, Rio Grande do Sul, Brazil; ⁵Oswaldo Cruz Foundation, Evandro Chagas Institute, Ananideua, Pará, Brazil; ⁶Tropical Medicine Foundation Dr. Heitor Vieira Dourado, Manaus, Amazonas, Brazil; ⁷Federal University of Rio de Janeiro, Institute of Chest Diseases, Clementino Fraga Filho University Hospital, Rio de

HIGHLIGHTS

- The Kit SIRE Nitratase® is the only commercial test using the Nitrate Reductase Assay
- It is a know-how transference from the public university to the national industry.
- It can be implanted in the health system in countries with high TB burden
- Highlighting the good accuracy, less time to results and less laborious.

Abstract: (1) Background: The Commercial Kit SIRE Nitratase® PlastLabor, is a drug susceptibility test kit used to detect *Mycobacterium tuberculosis* resistance to first-line TB treatment drugs. The present study aimed at evaluating its performance in a multicenter study. (2) Methods: To determine its accuracy, the proportion methods in Lowenstein Jensen medium or the BACTEC™MGIT™960 system was used as a gold standard. (3) Results: The study revealed that the respective accuracies of the kit with 190 *M. tuberculosis* clinical isolates, using the proportion methods in Lowenstein Jensen medium or BACTEC™MGIT™960 system as a gold standard, were 93.9% and 94.6%, 96.9% and 94.6%, 98.0% and 97.8%, and 98.0% and 98.9%, for streptomycin, isoniazid, rifampicin, and ethambutol, respectively. (4) Conclusion: Thus, the kit can rapidly screen resistance to streptomycin, isoniazid, rifampicin, and ethambutol. Additionally, it does not require sophisticated equipment; hence, it can be easily used in the laboratories of low and middle income countries.

Keywords: tuberculosis, MDR, nitrate reductase, technology

Table 1. Drug susceptibility testing using Kit SIRE Nitratase®; sensitivity, specificity, and accuracy using the Lowenstein Jensen medium as a standard (n = 98)

Drug	Drug susceptibility profile by NRA	Proportion Method using Lowenstein Jensen		Sensitivity (%)	Specificity (%)	Accuracy (%)
		Resistant	Susceptible			
Streptomycin ^a	Resistant	33	6	100	90.8	93.9
	Susceptible	0	59			
Isoniazid ^b	Resistant	41	3	100	94.7	96.9
	Susceptible	0	54			
Rifampicin ^c	Resistant	28	0	93.3	100	98
	Susceptible	2	68			
Ethambutol ^d	Resistant	17	0	89.5	100	98
	Susceptible	2	79			

Legend: a: Kappa = 0.87 (very good); b: Kappa = 0.94 (very good); c: Kappa = 0.95 (very good); d: Kappa = 0.93 (very good). Included sites: Pará/IEC, Manaus/FMT, IPEC/Fiocruz, IDT/UFRJ, and CDCT/RS.

Table 2. Drug susceptibility testing using Kit SIRE Nitratase®; sensitivity, specificity, and accuracy using the BACTEC™ MGIT™ 960 system as a standard (n = 92)

Drug	Drug susceptibility profile by NRA	BACTEC MGIT 960		Sensitivity (%)	Specificity (%)	Accuracy (%)
		Resistant	Susceptible			
Streptomycin ^a	Resistant	33	5	100	91.5	94.6
	Susceptible	0	54			
Isoniazid ^b	Resistant	38	0	88.4	100	94.6
	Susceptible	5	49			
Rifampicin ^c	Resistant	33	0	94.3	100	97.8
	Susceptible	2	57			
Ethambutol ^d	Resistant	25	1	100	98.5	98.9
	Susceptible	0	66			

Legend: *a*: Kappa = 0.89 (very good); *b*: Kappa = 0.89 (very good); *c*: Kappa = 0.95 (very good); *d*: Kappa = 0.97 (very good). Included sites: Adolfo Lutz/SP, CRPHF/Fiocruz, Dourados/UFGD, and IHMT/Lisboa.

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